

Certificate of Analysis AIHA-LAP EMLAP# 102977

43760 Trade Center Place Suite 100 Sterling, Virginia 20166 (877) 648-9150 www.aerobiology.net

Aerobiology Sample Report
43760 Trade Center Place
Sterling, Virginia 20166
Attn: John Doe
Project: USP 797 Counts Only Test

Date Collected: 06/14/2018
Date Received: 06/2018
Date Reported: 06/21/2018
Project ID: 18021276

Ornalition of Commission Descints Assemble

Condition of Sample(s) Upon Receipt: Acceptable

Client Sample #: 1 Lab Sample #: 18021276-001

Sample Location: Location 1

Test: 1204, USP 797 Culture, Surface Bacterial Counts Only

Results: 20 CFU/25cm² Area: 25 (cm²)

MRL: 1

Page 1 of 2

Client Sample #: 1 Lab Sample #: 18021276-001

Sample Location: Location 1

Test: 1206, USP 797 Culture, Surface Fungal Counts Only

Results: 13 CFU/25cm² Area: 25 (cm²)

MRL: **1**

Client Sample #: 2 Lab Sample #: 18021276-002

Sample Location: Location 2

Test: 1204, USP 797 Culture, Surface Bacterial Counts Only

Results: 18 CFU/25cm² Area: 25 (cm²)

MRL: **1**

Client Sample #: 2 Lab Sample #: 18021276-002

Sample Location: Location 2

Test: 1206, USP 797 Culture, Surface Fungal Counts Only

Results: 17 CFU/25cm² Area: 25 (cm²)

MRL: 1



Certificate of Analysis AIHA-LAP EMLAP# 102977

43760 Trade Center Place Suite 100 Sterling, Virginia 20166 (877) 648-9150 www.aerobiology.net

Aerobiology Sample Report
43760 Trade Center Place
Sterling, Virginia 20166
Attn: John Doe

Date Collected: 06/14/2018
Date Received: 06/15/2018
Date Analyzed: 06/20/2018
Date Reported: 06/21/2018
Project ID: 18021276

Project: USP 797 Counts Only Test

Condition of Sample(s) Upon Receipt: Acceptable

Project ID: 18021276

Page 2 of 2

USP 797 Class and Action Levels

ISO Clean Room Classification	ISO, 0.5 u/m ³ Particulate	Viable Air Sampling 400-1000 CFU/m ³	Surface Contact CFU/plate	Gloved Fingertip CFU/plate	Gloved Fingertip CFU/plate Gown Validation
Class 5	3,520	>1	>3	>3	>0
Class 7	352,000	>10	>5	N/A	N/A
Class 8 or Worse	3,520,000	>100	>100	N/A	N/A

Source PIC/S, 2007

Footnotes and Additional Report Information

- 1. Regardless of the number of CFU identified, further corrective actions are required if any pathogenic organisms are identified. It is therefore suggested to identify any colonies seen on the plate to genus level to rule out pathogens such as: gram-negative rods bacteria, and coagulase positive staphylococcus spp., yeasts, and mold.
- 2. Regardless of ISO Class, any fungal identification on an air or surface sample will cause the sample to be Out of Compliance.
- 3. Positive-hole correction factor is a statistical tool which calculates a probable count from the total raw count, taking into account multiple particles can impact on the same hole. For this reason the sum of calculated counts may be less than the positive hole corrected total.
- 4. TSA (Tryptic Soy Agar) for bacteria is incubated at 30-35°C for 2 days. MEA (Malt Extract Agar) or other suitable fungal media is incubated at 26 30°C for 5 to 7 days.
- 5. MEDIA CONTROLS. An unexposed TSA plate or MEA plate from each sampling event/project should be submitted for quality control purposes. The lot number for controls should be the same as those plates being submitted for analysis.
- 6. Semi-annual monitoring for viable bacteria and fungi in air, surface contact plates, gloved fingertip and particulates is required for both Class 5 and Class 7 defined areas.
- 7. Viable cultures must be collected using an impaction style sampler for volumetric capture. A sufficient volume of air (400 to 1000 liters) should be tested at each location to obtain the sensitivity and detection limit necessary for class action levels.
- 8. Standard contact plates have an area of 25 cm², unless otherwise noted in the sample area.
- 9. The results in this report are related to this project and these samples only.
- 10. **MRL** Units for USP 797 Cultures are as follows: AIR is CFU/m³, SURFACE is CFU/25cm², and CONTROL is colony/sample. **MRL**: Minimum Reporting Limit.
- 11. TARGET IDENTIFICATIONS: Any gram-negative rod, Staphylococcus aureus, yeast and molds
- 12. Non-sporulating colony is a colony of a filamentous mold on an agar plate that is not producing spores and/or conidiophores that allows the analyst to identify it further than a non sporulating colony. Identification structure must be present for identification.
- 13. If the final quantitative result is corrected for contamination based on the blank, the blank correction is stated in the sample comments section of the report.

Due to rounding totals may not equal 100%.

Suzanne S. Blevins, B.S., SM (ASCP) Laboratory Director

Syrace S. Bluing

© 2018 Aerobiology Laboratory Associates, Inc. All rights reserved.